

PULSE OXIMETER INSTRUCTION MANUAL

MD300C318T

Copyright

Our company owns all rights to this unpublished work and intends to maintain this work as confidential. Our company may also seek to maintain this work as an unpublished copyright. This publication is to be used solely for the purpose of reference, operation, maintenance or repair of our equipment. No part of this can be disseminated for other purposes.

In the event of inadvertent or deliberate publication, our company intends to enforce its right to this work under copyright laws as a published work. Those having access to this manual may not copy, disseminate or disclose the information in this work unless expressly authorized by our company.

All information contained in this manual is believed to be correct. Our company shall not be liable for errors contained herein nor for incidental or consequential damages in connection with the furnishing, performance, or use of this material. This publication may refer to information and protected by copyrights or patents and does not convey any license under the patent rights of our company, nor the rights of others. Our company does not assume any liability arising out of any infringements of patents or other rights of third parties.

Our Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Also, content of the manual is subject to change without prior notice.

© 2011 CHOICEMMED

All trademarks are the property of **CHOICEMMED** unless otherwise noted.

Bluetooth® and the Bluetooth® Logo ⌘ are registered trademarks of Bluetooth SIG, Inc.

V1.0318T

Issued Date 20-5-2011

Safety Information

Please read this manual thoroughly before using the pulse oximeter! Keep it in hand for future reference. Warnings alert the user to potential serious outcomes, such as injury or adverse events to the patient or user. Cautions alert the user to exercise care necessary for the safe and effective use of the pulse oximeter. Notes contain important information that may be overlooked or missed.

Warnings!

- DO NOT strike or needle the battery.
- Keep away from source of fire and/or heat.
- DO NOT disassemble the oximeter or its accessories.
- DO NOT use the pulse oximeter in an MRI or a CT environment.
- DO NOT use the pulse oximeter in the presence of flammable anesthetics.
- Explosion hazard: DO NOT use the pulse oximeter in an explosive atmosphere.
- Chemicals from a broken OLED panel are toxic when ingested. Use caution when the oximeter has a broken display screen.
- The pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
- Check the pulse oximeter application site frequently to determine the positioning of the measurement and circulation and skin sensitivity of the patient.
- Prolonged use or the patient's condition may require changing the measurement site periodically. Change measurement site and check skin integrity, circulatory status, and correct alignment at least every 4 hours.

- Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- Only the authorized service personnel can replace the battery or repair this device. This device uses a fixed lithium-ion battery inside. Do not try to replace the battery by yourself at any time. For longer battery life, only charge the battery when the battery power is empty.
- This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Cautions!

- Inaccurate measurements may be caused by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid may cause inaccurate readings.
- Significant levels of dysfunctional hemoglobin (such as carbonyl-hemoglobin or methemoglobin) may cause inaccurate readings.
- Intravascular dyes such as indocyanine green or methylene blue may cause inaccurate readings.
- SpO₂ measurements may be adversely affected in the presence of upper ambient light. Shield the sensor area (with a surgical tower, or direct sunlight, for example) if necessary.
- Excessive patient movement may cause inaccurate readings.
- Upper-frequency electrosurgical interference may cause inaccurate readings.
- Venous pulsations may cause inaccurate readings.
- Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line may cause inaccurate readings.
- The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia may cause inaccurate readings.

- Operation of the pulse oximeter may be affected by the use of an electrosurgical unit (ESU).
- The pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO₂ measurement. Verify that nothing is hindering the pulse measurement before replying on the SpO₂ measurement.
- The patient is in cardiac arrest or is in shock may cause inaccurate readings.
- Fingernail polish or false fingernails may cause inaccurate SpO₂ readings.
- Federal Law (U.S.A) restricts this device to sale by or on the order of a physician.
- All **CHOICEMMED** devices are designed to be compliant with rules and regulations in locations they are sold and will be labelled as required.
- Any changes or modifications to **CHOICEMMED** equipment, not expressly approved by **CHOICEMMED**, could void the user's authority to operate the equipment.

Contents

CHAPTER 1 General Information	7
1.1 Measuring Principle	7
1.2 Product Features	8
1.3 Intended Use	8
1.4 Appearance Introduction	8
1.5 Description of Symbols	10
CHAPTER 2 Settings	11
2.1 Configuring Your Pulse Oximeter	11
2.2 Selecting a Display Mode	11
CHAPTER 3 Take a Measurement	12
CHAPTER 4 Battery Charge	13
4.1 Power Supply	13
4.2 Battery Charge	13
CHAPTER 5 Data Transmission	14
CHAPTER 6 Specifications	16
CHAPTER 7 Maintenance and Calibration	18
CHAPTER 8 Troubleshooting	19

CHAPTER 1 General Information

Oxygen saturation is a percentage of oxyhemoglobin (HbO_2) capacity, compounded with oxygen, by all combinative hemoglobin (Hb) capacity in blood. In other words, it is consistency of oxyhemoglobin in blood. It is a very important parameter for the Respiratory Circulation System. Many respiratory diseases can result in low oxygen saturation in human blood. Additionally, the following factors can reduce oxygen saturation: Automatic regulation of organ dysfunction caused by Anesthesia, Intensive Postoperative Trauma, injuries caused by some medical examinations. That situation might result in light-headedness, asthenia, and vomiting. Therefore, it is very important to know the oxygen saturation of a patient so that doctors can find problems in a timely manner.

The pulse oximeter integrated with Bluetooth® wireless technology is designed to monitor vital signs, such as SpO_2 and pulse rate anytime, anywhere.

1.1 Measuring Principle

An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin(HbR) and Oxyhemoglobin (HbO_2) in red light and near-infrared light zones. The photoelectric oxyhemoglobin inspection technology is adopted in accordance with capacity pulse scanning and recording technology, so that two beams of different wavelength of lights (660nm red light and 940nm near infrared light) can be focused onto human nail tip through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on the display through process in electronic circuits and microprocessor.

Diagram of Operation Principle:

1. Red and Infrared light emission diodes
2. Red and Infrared photodiode



1.2 Product Features

- ▶ Color OLED screen displays SpO₂, PR, waveform and pulse bar; 4 display modes.
- ▶ Rechargeable lithium-ion battery, battery charging via USB cable or adapter.
- ▶ Real-time and non real-time Bluetooth® data transmission.
- ▶ Real-time inner clock.
- ▶ Software upgrade with Bootloader.
- ▶ Suitable for pediatrics and adults.

1.3 Intended Use

The MD300C318T Fingertip Pulse Oximeter is intended for continuous use or spot checking in measuring and displaying functional arterial oxygen saturation (SpO₂) and pulse rate of patients in hospitals and home care. It is intended for adult and pediatric patients on finger between 0.3-1.0 inch (0.8 - 2.5 cm) thick.

1.4 Appearance Introduction




Fig.1.2



Fig.1.3

Description of the Device

1. Display screen: provides an organic light-emitting diode (OLED) display panel. [Figure 1.2]
2. Power switch & function button is used to turn on the device and set the display orientation [Figure 1.2-2]. See full description of the functions below, under items 1-10.
3. USB cable port: is used to charge the device and transfer data. [Figure 1.3]

The display screen is illustrated in the figure below




Fig1.4













Description of fig1.4

- 1) Time display: The time shown in the figure is 11:45.
- 2) SpO₂% refers to the oxygen saturation in blood.
- 3) Tabulation of the measured SpO₂% reading. Here, it is valued at 98%.
- 4) This field indicates SpO₂% Plethysmograph signal's dimensions.
- 5) This field indicates the measurement of pulse rate (PR). Here, the PR is valued at 60 beats per minute (bpm).

- 6) The pulse graph bar indicates pulse amplitude.
- 7) PRbpm refers to Pulse Rate beats per minute.
- 8) This field indicates the level of battery power in the pulse oximeter. When battery power is low, the icon will be empty and turn to red.
- 9) Bluetooth® transmission indicator.
- 10) Brightness level indicator.

NOTE: The illustrations used in this manual may differ slightly from the appearance of the actual product.

1.5 Description of Symbols

Symbols	Descriptions	Symbols	Descriptions
	Power or function button		Bluetooth® Transmission Indicator
	Power low indicator (red)		Battery power indicator
	Non AP equipment		European Union Approval
	Type BF Applied Part		Resistant to liquid ingress
	Date of Manufacture		Manufacturer
SN	Serial Number		Attention! See Instructions for Use
	Authorised representative in the European community		

CHAPTER 2 Settings

2.1 Configuring Your Pulse Oximeter

There are some settings which can be adjusted by software, including date and time, switching, real-time and non real-time data transmission modes, powering off, and updating the software of the pulse oximeter.

2.2 Selecting a Display Mode

You can adjust the display mode of your oximeter when performing SpO₂ and pulse rate measurements, every time you press the button, the pulse oximeter will switch to another display mode.

During the measuring, each time you press the button the oximeter will switch to another display mode, there are 4 display modes shown as below:




Fig 2.1

CHAPTER 3 Take a Measurement

1. Open the clamp and insert a finger into the oximeter, as illustrated in the right figure. Then gently release the clamp. Make sure that the distance between you and other Bluetooth® device within 10m. (30 ft.)




Fig3.1

2. Press the power button, the display will light up and after a few seconds, the device will display SpO₂% and pulse rate measuring values. The oximeter will also transmit the recorded information via Bluetooth® to other Bluetooth® device. After transmitting data to other Bluetooth® device, the oximeter will shut down automatically.

NOTES:

- ▶ Keep your tested hand still during the measurement.
- ▶ Do not pour or spray liquids onto the oximeter, and do not allow any liquid to enter any openings in the device. Make the oximeter dry thoroughly before reusing.
- ▶ The oximeter will shut down after a period of time if no measurement is being taken.
- ▶ Clean the rubber touching the finger probe inside of the pulse oximeter with a soft cloth dampened with 70% isopropyl alcohol, and cleans the test finger probe using alcohol before and after each measurement.

CHAPTER 4 Battery Charge

4.1 Power Supply

Battery Model: SP080

Battery Type: One 3.7V Lithium Ion Rechargeable Battery

4.2 Battery Charge

4.2.1 Connect the oximeter with the attached USB cable, shown as the following figure.

4.2.2 Connect the other end of the USB cable to the attached charger and plug the charger into an electrical outlet. Charging will complete in 4-5 hours. The battery indicator will stop flickering if the charging process is completed.




Fig.4.1 Connection




Fig.4.2 Charging finished

Charging temperature: 0°C~40°C

Charging Voltage: 110Vac~240Vac

 **Warning!**

- ▶ Avoid strongly impacting to the oximeter.
- ▶ Avoid the oximeter to be exposed straightly to sunlight.
- ▶ Only use the attached battery charger. Otherwise, damage even danger may be caused to oximeter or person.
- ▶ Do not charge the battery for more than 5 hours, or that may cause damage to the battery.
- ▶ Do not charge or preserve the battery in too hot or cold environment, the proper temperature is 0~40℃ for charging. Please refer to details in **Chapter 6 Specifications**.
- ▶ Charging is not recommended when using the oximeter.

CHAPTER 5 Data Transmission

The measurement results saved in the pulse oximeter can be uploaded to other Bluetooth® devices for review or management.

The modes of data transmission include Realtime Bluetooth® transmission and Non-realtime Bluetooth® transmission. Please refer to "MD300C318T Communication Protocol" for details of data transmission.

• Realtime Bluetooth® transmission(default):

- ▶ Bluetooth® icon ✕ will flash when MD300C318T is not connected.
- ▶ Bluetooth® icon ✕ will not flash when the valid pairing of MD300C318T is connected.
- ▶ If Bluetooth® is not connected, the measured data will not be saved and MD300C318T will be in a wait state .(If Bluetooth® has not been connected for 3 minutes, MD300C318T will shutdown automatically.)
- ▶ If Bluetooth® connection is successful, MD300C318T will upload the values of blood oxygen saturation, pulse rate, blood oxygen volume wave, battery power to other Bluetooth® devices; Meanwhile, MD300C318T can respond and execute commands from other Bluetooth® devices. (Please see MD300C318T Communication Protocol.)

• Non-realtime Bluetooth® transmission:

The data transmission screen is shown in Fig.5.1.




Fig.5.1

- ▶ MD300C318T will automatically judge the effective measured data (6 single point effective measurements will be stored).
- ▶ Bluetooth® will automatically open and upload the measured data after measurement.
- ▶ MD300C318T software system can set the number of stored data of one successful measurement, the number ranging from 1-10.
- ▶ If data upload fails, the measured data will be saved automatically and uploaded with next measured data. At most 100 pieces of failed-uploaded data can be stored and re-uploaded.
- ▶ The distance between the pulse oximeter and other Bluetooth® device should be not more than 10m.
- ▶ The measurement data is saved every 4s and at most 100 pieces of data records can be saved in the pulse oximeter.

CAUTION: The Bluetooth® transmission may be influenced in different environments, such as walls, metallic doors, steel wire netting and an MRI or a CT environment and so on.

CHAPTER 6 Specifications

Display:

Display type: OLED display
SpO₂ display range: 0~99%
SpO₂ measuring range: 70~99%
PR display range: 0~235 BPM
PR measuring range: 30~235 BPM
PR display mode: Amplitude Bar
Data update period: <15 s

Resolution:

SpO₂%: 1%
Pulse rate: 1BPM

Measurement Accuracy:

SpO₂: 70%~99% ±3%; ≤69% unspecified.
PR: 30~99 bpm ±2 bpm; 100~235 bpm ±2%.

Antenna Information:

Antenna Type/Patten: Internal
Frequency Range: 2402 to 2480 MHz

LED Wavelengths and Output Power:

Red: approximately 660nm @0.8mW maximum average
Infrared: approximately 940nm @0.8mW maximum average

Bluetooth®

The device name format of MD300C318T is "B *****". ("*****" denotes numbers from 0 to 9, for example B110500001.) The device name will be labeled on the back of MD300C318T pulse oximeter.

Wireless transmission distance: 0~10m
Bluetooth® default paired password: 0000.

Physical Characteristics:

Dimensions: 60mm x 36mm x 38mm (2.36"x 1.42"x 1.50")
Weight: 50g (including the lithium ion battery)

Battery Life:

Typical lithium ion battery: 300 cycles
Charge temperature: 0°C~40°C

Environment Conditions:

Operating Temperature: 5°C~40°C
Storage Temperature: -20°C~-55°C
Humidity: ≤80%, no condensation in operation
≤93%, no condensation in storage

Accessories:

AC/DC Adapter	1 piece
Instruction Manual	1 piece
Lanyard	1 piece
USB cable	1 piece

CHAPTER 7 Maintenance and Calibration

- 7.1 Please charge the battery in time when the low power indicator flickering.
- 7.2 Clean the surface of the fingertip oximeter before it is used in diagnosis for patients.
- 7.3 It is better to preserve the product in the environment specified in this manual.
- 7.4 It is recommended that the product should be kept in a dry environment. A wet ambient might affect the lifetime, even damage the device.
- 7.5 Please follow the law of local government to deal with the used oximeter and the accessories.

CHAPTER 8 Troubleshooting

Problems	Reasons	Solutions
SpO ₂ % or pulse rate does not display.	<ol style="list-style-type: none"> 1. Finger is not plugged correctly. 2. Patient's SpO₂ value is too low to be measured. 	<ol style="list-style-type: none"> 1. Retry by plugging the finger. 2. There is excessive illumination. 3. Measure other patients to make sure that no problem exists in the product. Go to a hospital in a timely manner for an exact diagnosis.
SpO ₂ % or PR is shown unstably.	<ol style="list-style-type: none"> 1. Finger might not be plugged deep enough into the clamp probe. 2. Excessive patient movement. 	<ol style="list-style-type: none"> 1. Retry by inserting the finger to the end. 2. Stop moving the finger, hand or body.
The oximeter cannot be powered on.	<ol style="list-style-type: none"> 1. No battery or low power of battery. 2. Battery might be installed incorrectly. 3. The Monitor might be damaged. 	<ol style="list-style-type: none"> 1. Please replace battery. 2. Please reinstall the battery. 3. Please contact with local customer service centre.
Display suddenly turns off.	<ol style="list-style-type: none"> 1. The oximeter is automatically powered off when no signal is detected longer than 8 seconds. 2. The batteries power is too low to work. 	<ol style="list-style-type: none"> 1. Relocate the probe on another finger or restart the oximeter and be sure the signal strength is strong enough for stable display. 2. Replace the battery.
Error 1 Error 2 Error10	ROM error RAM error EEPROM damaged or dry joint	Please contact with local customer service centre. Please contact with local customer service centre. Please contact with local customer service centre.

CHAPTER 9 Declaration

FCC-ID: WWIMD300C318T

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference which will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Ask the dealer or an experienced radio/TV technician for help.

EMC of this product comply with IEC60601-1-2 standard

Guidance and Manufacture's declaration – electromagnetic emissions-For all EQUIPMENT and SYSTEMS

Guidance and Manufacture's declaration - electromagnetic emission

The MD300C318T PULSE OXIMETER is intended for use in the electromagnetic environment specified below. The customer of the user of the MD300C318T PULSE OXIMETER should assure that it is used in such and environment.

Emission test	Compliance	Electromagnetic Environment – guidance
RF emissions CISPR 11	Group 1	The MD300C318T PULSE OXIMETER uses RF energy only for internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The MD300C318T PULSE OXIMETER is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and Manufacture's declaration – electromagnetic immunity- For all EQUIPMENT and SYSTEMS**Guidance and Manufacture's declaration - electromagnetic immunity**

The MD300C318T PULSE OXIMETER is intended for use in the electromagnetic environment specified below. The customer of the user of the MD300C318T PULSE OXIMETER should assure that it is used in such and environment.


Immunity test	IEC 60601 test level	Compliance Level	Electromagnetic Environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	+/- 6kV contact +/- 8kV air	+/- 6kV contact +/- 8kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	+/-2Kv for power supply lines	+/-2Kv for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/-1Kv differential mode	+/-1Kv differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles 5% U_T (>95% dip in U_T) for 5 sec.	<5% U_T (>95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles 5% U_T (>95% dip in U_T) for 5 sec.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the MD300C318T PULSE OXIMETER requires continued operation during power mains interruptions. It is recommended that the MD300C318T PULSE OXIMETER be powered from an uninterruptible power supply.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level

Guidance and Manufacture's declaration – electromagnetic immunity- For EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and Manufacture's declaration - electromagnetic immunity

The MD300C318T PULSE OXIMETER is intended for use in the electromagnetic environment specified below. The customer of the user of the MD300C318T PULSE OXIMETER should assure that it is used in such and environment.

Immunity test	IEC 60601 test level	Compliance Level	Electromagnetic Environment – guidance
Radiated RF IEC 61000-4-3	3 V/m 150 KHz to 80 MHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the MD300C318T PULSE OXIMETER, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{V_t} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_t} \right] \sqrt{P} \quad 80\text{MHz to } 800 \text{ MHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \left[\frac{7}{E_t} \right] \sqrt{P} \quad 800\text{MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with following symbol: </p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations, Electromagnetic propagation is affected by absorption and reflection structures, objects and people.

Field strengths from fixed transmitters, such as base station for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MD300C318T PULSE OXIMETER is used exceeds the applicable RF compliance level above, the MD300C318T PULSE OXIMETER should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the MD300C318T PULSE OXIMETER. Over the frequency range 150kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEMS - For all EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the MD300C318T PULSE OXIMETER

The MD300C318T PULSE OXIMETER is intended for use in electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MD300C318T PULSE OXIMETER can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MD300C318T PULSE OXIMETER as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	80 kHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 kHz to 2.5 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	3.69
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

 Beijing Choice Electronic Technology Co.,Ltd.

BailangyuanB1127-1128,Fuxing R,A36

100039 Beijing

PEOPLE'S REPUBLIC OF CHINA

EC	REP
----	-----

 Eiffestraße 80,20537
Hamburg GERMANY

Shanghai International
Holding Corp. GmbH (Europe)

